



Terry Clothing Protector

The Facts:

- Crafted from a warp knit terry blend of 90% cotton and 10% polyester
- Soft and absorbent surface for maximum comfort
- Available with a 100% impermeable white barrier, ensuring no liquids pass through to protect clothing effectively
- Edged with resilient 65/35 cotton/polyester woven binding for enhanced durability and longevity



Contact your sales consultant



Also available on the MIP store
store.mipinc.com



Specifications

IMPERMEABLE

STYLE#	DIMENSIONS	WEIGHT	COLOUR
30836/SB	18" x 36"	7.3 oz	Blue
30836/SW	18" x 36"	7.3 oz	White

WITH SNAPS

STYLE#	DIMENSIONS	WEIGHT	COLOUR
40836/SB	18" x 36"	4.7 oz	Blue
40836/SW	18" x 36"	4.7 oz	White

MIP

The Softer Side of Healthcare