



Microfiber Isolation Gown


2021 v. 2


The Facts:





- Reversible style and stain resistant properties to reduce labour time and costs
- Durable and quick-dry including knit cuffs
- Fabric meets level 2 requirements
- Carbon yarns reduce static build-up



STYLE#	EHX4200/WA	GOWN/PPE-3XL (Back tie/Non-reversible)
FULL SWEEP	63"	74"
SLEEVE INSEAM	24"	24"
WEIGHT	280 gr	307 gr

 Contact your sales consultant

 Also available on the MIP store
store.mipinc.info

MIP

The Softer Side of Healthcare