



# Cotton Bed Pads

## The Facts:

- All these pads feature a blend of cotton and polyester for the surface material, providing a comfortable and durable surface, with a polyester/rayon soaker for high absorbency and a Vintex barrier for effective moisture protection
- The quilt pattern enhances both the aesthetic and functional design, offering a gentle and soothing texture while managing moisture efficiently and protecting against spills and stains
- Available in a range of sizes, from standard to larger dimensions, to suit different bed types



Contact your sales consultant



Also available on the MIP store  
[store.mipinc.com](http://store.mipinc.com)



# MIP

The Softer Side of Healthcare



# Specifications

## WHOOTEX BED PAD

STYLE #	WEIGHT	DIMENSIONS	COLOUR
526/SQ	34"x 36"	429 g	White
527/Q	50"x 60"	1205 g	White
528/Q	34"x 45"	587 g	White

## COTTON BED PAD

STYLE #	WEIGHT	DIMENSIONS	COLOUR
500/Q	34"x 36"	508 g	White
508/Q	34"x 45"	677 g	White

## COTTON CHAIR PAD

STYLE #	WEIGHT	DIMENSIONS	COLOUR
515/Q	17"x 22"	172 g	White



# MIP

The Softer Side of Healthcare