



# FUSION UNDERPADS

## Patient care underpads

### The Facts


- “Stay Dry” technology features distinct characteristics to promote skin care: rapid wicking and high absorbency
- Ultra thin profile lays flat
- Reduces pressure points and eliminates wrinkling & surface puckering
- Ensures fewer bed linen changes
- An eco-friendly alternative to disposable incontinence pads



FUSION creates a smooth surface for a care recipient's sensitive skin. With skin integrity at the forefront of its mandate, FUSION helps protect the patient's skin from moisture, shear and friction.

## Specifications

STYLE#	FORMAT	DIMENSION	WEIGHT	SURFACE COLOR	BARRIER
9100/FP	Standard	35" x 35"	16.4 oz	Turquoise	Blue PU
9118/FP	Bariatric	34" x 45"	19.4 oz	Turquoise	Blue PU
9105/FP	Chairpad	17" x 24"	5.4 oz	Turquoise	Blue PU



Contact your sales consultant



Also available on the MIP store  
[store.mipinc.info](http://store.mipinc.info)



# MIP

The Softer Side of Healthcare

2023 v.1

HEAD OFFICE: 9300 Rue du Parcours, Montreal (Quebec) H1J 2Z1 Canada • Tel.: (514) 356-1224 | Toll Free: 1 (800) 361-4964 | Fax: (514) 356-0055

[mip.ca](http://mip.ca)