

# INTEGRIDERM<sup>®</sup>

PRESSURE REDISTRIBUTION MATTRESSES  
BY MIP

**Putting Integrity & Comfort First**




**MIP**   
[www.mipinc.com](http://www.mipinc.com)

**RESPURTEX® 4-WAY STRETCH**



Breathable, moisture wicking material ensures a healthy micro-climate for the skin while protecting against shearing

**MAXIMIZED SLEEPING SURFACE**



Enhances comfort when sleeping on or near the edge of mattress leading to better rest

**CROSS DIE-CUT**



Zoned cut-outs to enhance pressure redistribution and care recipient comfort

**STANDARD FIRE SHIELD**



Safety first. Fire-resistant sleeve around foam components to protect against flame

**HEEL SLOPE SUPPORT**



Gradual downward slope to decrease load on the heels to prevent skin breakdown

# The Integriderm Story

PRODUCT INTEGRITY, SKIN INTEGRITY, CARE INTEGRITY



**CARE RECIPIENTS**

- Enhanced comfort for better rest
- Provides a happy and healthy skin environment



**STAFF**

- Promotes staff safety with our innovative design (handles out of sight, low slip bottom fabric)
- Minimal maintenance required
- Healthy skin and comfortable patients help staff prioritize their care efforts



**BUDGET & ADMINISTRATION**

- Affordable quality
- Offers peace of mind (avoid replacement costs with warranties up to 7 years)

MIP Inc. is dedicated to offering an economic line of high quality mattresses to meet the needs of a care recipient for all periods of time. Our Integriderm formula reduces the occurrence of pressure injuries while providing better quality rest through an exceptional design.


## Pressure Injuries – The Whole Story

**WHAT ARE THEY?**


Pressure injuries result from extended pressure interrupting blood flow to the skin, causing problems ranging from discoloration, loss of sensation and scarring to full-fledged skin loss and soft tissue damage.

**HOW DO THEY OCCUR?**

They can result from pressure alone or the combination of pressure and friction or shearing.




**WHO IS AT RISK?**



**HOW CAN THEY BE PREVENTED?**

Pressure injuries are easier to prevent than to heal:

- Using pressure redistribution devices to alleviate pressure from key points
- Checking skin daily for discoloration
- In most cases, ensuring the head of the bed is raised no more than 30° and for short periods, if possible



# Integriderm Comparative Mattress Features

	MIP9000V/BR	MIP9080	MIP6080	MIP8080
 CROSS DIE-CUT VISCO-ELASTIC GEL INFUSED TOPPER	+			
 FOAM LAYERS	4	3	2	3
 RESPURTEX 4-WAY STRETCH COVER	+	+	+	+
 MAXIMIZED SLEEPING SURFACE	+	+	+	+
 CROSS DIE-CUT DESIGN	+	+	+	+
 STANDARD FIRE SHIELD	+	+	+	+
 HEEL SLOPE SUPPORT	+	+	+	+
 FULLY CLOSED COVER	+	+	+	+
 ANTI-SKID BASE	+	+	+	+
 BOLSTERED RAILS	+	+	+	+
 TECHNOSTAT ANTI-BACTERIAL FILTER	+	+	+	+
 MOISTURE RESISTANT WELDED HANDLES	+	+		+
 WARRANTY				



## Integriderm:

### AN INTEGRAL ELEMENT OF THE REVOLUTION BED BY MIP

A properly made bed enhances the function of pressure redistribution mattresses.

The Revolution Bed concept was designed with the care recipient and the care giver in mind. Take the first step with an Integriderm mattress and learn more about the *Revolution Bed* make-up today.

### THE REVOLUTION BED



### Putting Integrity & Comfort First



CALL US  
1 (800) 361-4964



EMAIL US  
info@mipinc.com



GO TO OUR WEBSITE  
mipinc.com